



CHECKMATE 2™

Simple and Reliable Daily QA

A therapy beam constancy meter intended for simple radiation output checks.

The CHECKMATE 2 Advantage

- ▶ Simple to use
- ▶ Cost effective solution
- ▶ Automatic temperature and pressure correction
- ▶ Accurate measurement of:
 - > CAX output
 - > Repeatability
 - > Monitor chamber linearity
- ▶ Built-in rechargeable batteries
 - > Safely keep cables off the treatment floor

Hardware

- ▶ Single vented ionization chamber at a depth of 1.0g/cm²
- ▶ Precision alignment template with light field markings at 10 x 10cm and 20 x 20cm
- ▶ Large, easy to read display
 - > Easily viewed at the control station bunker monitor
- ▶ User calibration
 - > Stores up to 15 calibration values
 - > Automatic memory position advance


Process

- 1 STEP** Turn on the power switch.
- 2 STEP** Place CHECKMATE 2 on the treatment couch and align the crosshairs and the light field to the instrument template.
- 3 STEP** Expose CHECKMATE 2 to the beam.
- 4 STEP** Read the measured dose on the numeric display on the end of the instrument.



The numeric dose value displayed indicates the percentage of the stored calibration value. The displayed value can be viewed on the bunker monitor or directly by entering the vault.

Features and Specifications



Detector type:	Vented Ion Chamber	Short-term reproducibility:	0.2%
Detector quantity:	1 - CAX	Long-term reproducibility:	0.5%/year
Detector volume (cm ³):	0.6	Temperature / Pressure correction:	Automatic
Field size (cm):	20.0 x 20.0	Calibration factor memory positions:	15
Inherent buildup (g/cm ²):	1.0 ± 0.1	Display units:	Percent of calibration value; set by user
Inherent backscatter (cm):	Acrylic, 3.6	Display range:	0.1% to 999.9%
Radiation measured:	<ul style="list-style-type: none"> • Electrons, 6MeV to 25MeV • Photons, Co-60 to 25MV 	Dimensions / Weight:	25.5 x 40.0 x 5.5cm/4.9kg
Beam limits:	<ul style="list-style-type: none"> • Max. dose per pulse: 3.65cGy • Max. avg. dose rate: 3410cGy/min 	Battery operation:	8 hours minimum
		Number of connection cables:	Power charger for internal batteries

All data used is best available at time of publication. Data is subject to change without notice. All Content ©2010, Sun Nuclear Corporation. All Rights Reserved.