

# Sun Nuclear: Beam quality correction factors for radiation therapy with high energy photon beams



Ionization chamber measurements of the absorbed dose to water in high photon energy beams are an essential part of radiation therapy dosimetry. The response deviation of an ionisation chamber in reference conditions and clinical beam qualities is corrected for by the beam quality correction factor  $k_Q$ , which provides one of the most significant sources of uncertainty in the measurement of the absorbed dose to water.

$k_Q$  values depend on the design, size and material composition of the components and are chamber specific. International dosimetry protocols often provide tabulated data for  $k_Q$  values. However, the recent release of ICRU 90 [1] provided updated values of the mean ionization energy of water and graphite (I-values), as well as updated density correction parameters. Consequently,  $k_Q$  values based on ICRU 37 [2] should be revised.

While Monte Carlo simulations using science grade algorithms like EGSnrc [3] (National Research Council, Ottawa, Canada) are the method of choice to calculate kQ values for clinically used ionization chambers, the implementation and validation of ionization chamber models is not trivial. To ensure the highest quality standard, we connected two leading research groups, the Canadian national metrology institute NRC, and the Institute of Medical Physics and Radiation Protection of the University of Applied Sciences Mittelhessen, Germany, for the revision and calculation of  $k_Q$  according to ICRU 90 [1]. Both groups independently created and validated Monte Carlo models of the reference class ionization chambers SNC 125c and SNC 600c to determine  $k_Q$  values according to TRS 398 [4], AAPM TG 51 [5] and DIN 6800-2 [6]. The Monte Carlo calculated data sets agreed with experimentally determined kQ values presented in NRC reports PIRS-3224 [7] and PIRS-3327 [8].

Figure 1 presents Monte Carlo calculated beam quality correction factors kQ for the reference class ionization chambers SNC 125c and SNC 600c, independently determined by NRC and THM according to the TRS 398 [4] protocol. The kQ values were calculated for various clinical photon beams with nominal energies ranging from 4 MV up to 24 MV. A polynomial function was fitted to the joint data sets calculated by THM as well as NRC.




Figure 1: Calculated beam quality correction factors for the reference class ionization chambers SNC 125c and SNC 600c. Experimental data according to TRS 398 [4] is shown for the SNC 600c and available for the SNC 125c according to TG 51 [5] in Alissa et al. (2021) [9].

$k_0$  data according to TRS 398 [4], AAPM TG 51 [5] and DIN 6800-2 [6] for Sun Nuclear's ionization chambers SNC 125c and SNC 600c, as well as the corresponding fit parameters can be openly accessed in Alissa et al. (2021) [9].

## References

- [1] Seltzer S, Fernandez-Varea J, Andreo P, et al. Key data for ionizing-radiation dosimetry: measurement standards and applications, *ICRU Report 90*. 2016;
- [2] Berger M, Inokuti M, Anderson H, et al. Report 37. Journal of the International Commission on Radiation Units and Measurements. 1984;(2)
- [3] Kawrakow I, Rogers D, Mainegra-Hing E, Tessier F, Townson R, Walters B. EGSnrc toolkit for Monte Carlo simulation of ionizing radiation transport. 2000; doi:[10.4224/40001303](https://doi.org/10.4224/40001303)
- [4] IAEA. *TRS 398: Absorbed Dose Determination in External Beam Radiotherapy*. 2006.
- [5] McEwen M, DeWerd L, Ibbott G, et al. TG 51 Addendum to the AAPM's TG-51 protocol for clinical reference dosimetry of high-energy photon beams. *Med Phys*. Apr 2014;41(4):041501. doi:[10.1118/1.4866223](https://doi.org/10.1118/1.4866223)
- [6] DIN. Dosismessverfahren nach der Sondenmethode für Photonen- und Elektronenstrahlung Teil 2: Dosimetrie hochenergetischer Photonen- und Elektronenstrahlung mit Ionisationskammern. *DIN 6800-2*. Berlin: Beuth; 2020.
- [7] Muir B, McEwen M. Characterization of SNC 125c™ ionization chambers for measurements in linear accelerator photon beams. 2021. PIRS-3224.
- [8] McEwen M, Muir B. Characterization of SNC 600c™ ionization chambers for measurements in linear accelerator photon beams. 2021. PIRS-3327.
- [9] Alissa M, Zink K, Tessier F, Schoenfeld AA, Czarnecki D. Monte Carlo calculated beam quality correction factors for two cylindrical ionization chambers in photon beams. *Phys Med*. Dec 28 2021;94:17-23. doi:[10.1016/j.ejmp.2021.12.012](https://doi.org/10.1016/j.ejmp.2021.12.012)



**Mohamad Alissa, M.Sc.**

Mohamad Alissa is a clinical medical physicist at the University Hospital of Giessen and Marburg, Germany, and a PhD candidate at the Institute of Medical physics and radiation protection (IMPS) of the University of Applied Sciences Mittelhessen, Germany. He holds a master's degree in Radiation Protection and specializes in the field of radiation dosimetry in the presence of magnetic fields.



**Dr. Damian Czarnecki**

Damian Czarnecki is a post-doctoral research scientist at the Institute of Medical Physics and Radiation Protection (IMPS) of the University of Applied Sciences Mittelhessen, Germany. Damian holds a doctorate degree in physics and specialises in the fields of photon dosimetry, detector physics and Monte Carlo simulations.



**Dr. Andreas Schönfeld**

Andreas Schönfeld is a research associate at Sun Nuclear, the leader in quality assurance solutions for radiation therapy and diagnostic imaging. Andreas is based in Germany and supports product research, supervises the academic outreach program, and represents Sun Nuclear in DIN and IEC standardisation committees. He holds a doctorate degree in physics and is a certified clinical medical physicist.