




v6.0 Feature Set

With version 6.0 of the SunCHECK® Platform, we are pleased to introduce key updates, including Daily QA™ 4 Pro integration and SunCHECK Local.

See below for an overview of key v6.0 features and their areas of impact.



Feature	Workflow	Reporting	Analytics	Innovation	User-Driven
NEW! SunCHECK Local: Flexibility to install and use SunCHECK directly on a laptop or workstation, or enterprise-wide, to fit department and IT requirements	X				X
Mongo Database Removal: SunCHECK has moved from MongoDB to SQL for stronger data consistency and integrity.					X
Security & IT Improvements: Supports enterprise authentication and improved user access control					X
NEW! Daily QA™ 4 Pro: Enhanced daily machine QA in a single, indexed device Provides comprehensive QA coverage within a single setup, including output and profile constancy, photon/electron energy, field size/shift and imaging/positioning.	X	X		X	X
Trajectory Log Support: Supports Varian Medical Systems Ethos™ and Halcyon™ System Trajectory Log files, version 5.0		X	X		X

Platform-level features

Patient QA-focused features

Machine QA-focused features



Learn more about SunCHECK v6.0 on the **Release Highlights** page at sunnuclear.com



sunnuclear.com